



Calibration Kits for HF Instrumentation

BENCH TOP METERS

MICRO 100

#8219 White Light

Kit includes 0.0, 10 and 1000 NTU standards

- SKU #85529 – Primary Kit
- SKU #85555 – Deluxe Kit

MICRO 1000

#8268 White Light

Kit includes 0.0, 10 and 1750 NTU standards

- SKU #85431 – Primary Kit
- SKU #86186 – Deluxe Kit

#8221 Infrared

Kit includes 0.0, 10, 100 and 1750 FNU standards

- SKU #85531 – Primary Kit
- SKU #85557 – Deluxe Kit

ONLINE METER

MICRO TOL 3

#8276

Kit includes 0.0, 10 and 100 NTU standards

- SKU #86327 – Primary Kit
- SKU #86328 – Deluxe Kit

PORTABLE METERS

MICRO TPI

#8275

Kit includes 0.0, 10, 100 and 1000 NTU standards

- SKU #86325 – Primary Kit
- SKU #86326 – Deluxe Kit

MICRO TPW

#8252

Kit includes 0.0, 10, 100 and 1000 NTU standards

- SKU #86220 – Primary Kit
- SKU #86226 – Deluxe Kit

DRT 15CE

#8265

Kit includes 0.0, 10, 100 and 1000 NTU standards

- SKU #86169 – Primary Kit
- SKU #86170 – Deluxe Kit

PRODUCTS – PROVEN & TIME TESTED

- Non-toxic; safe handling during use
- Traceable to NIST Primary Standards
- Lot to lot accuracy of 1%
- Easy-to-use: No mixing or dilutions
- Shelf life guaranteed a year from DOS*
- Custom values available upon request



*DOS: Date of Shipment



Individual Standards for HF Instrumentation

We offer Individual Primary Standards in: 125mL, 250mL, 500mL, 1L, and 1 Gallon.
We also Offer: Sealed Secondary Standards.

HF White Light Instruments


Parent Bulk Item #	NTU
8000	0
8001	0.1
8002	0.2
8003	0.4
8004	0.5
8005	0.6
8006	0.8
8007	1.0
8008	2.0
8009	3.0
8010	4.0
8011	5.0
8012	6.0
8013	8.0
8014	10
8015	20
8016	30
8017	40
8018	50
8019	60
8020	80
8021	100
8022	200
8024	300
8025	400
8026	500
8027	600
8028	800
8029	900
8030	1000
8079	1750
8031	2000
8358	10,000

HF Infrared Instruments

Parent Bulk Item #	FNU
8000	0
8076	10
8329	100
8078	1000
8359	1750
8532	10,000

GFS CHEMICALS ALSO OFFERS:

- pH Buffers
- Conductivity Standards
- AA/ICP Standards
- UV/Vis Spectrophotometer Standards
- High Purity Acids
- High Purity Solvents



Updated 1/12/24