



Calibration Kits for Hach Instrumentation

BENCH METERS

2100 N/AN (TL2300/TL2350)

#8215

Kit includes 0.0, 20, 200, 1000, 4000 NTU standards

- SKU #85525 – Primary Kit
- SKU #85405 – Secondary Kit
- SKU #85551 – Deluxe Kit

2100 N IS/AN IS (TL2310/TL2360)

#8195

Kit includes 0.0, 20, 200, 1000, 4000 NTU standards

- SKU #86096 – Primary Kit
- SKU #86517 – Secondary Kit
- SKU #86097 – Deluxe Kit

TU5200 EPA

* Individual Standards Available

LEGACY METERS

2100 A

#8202

Kit includes 0.1, 0.5, 5.0 NTU standards

- SKU #85512 – Primary Kit
- SKU #85536 – Deluxe Kit

RATIO INSTRUMENTS

#8212

Kit includes 1.8, 18 and 180 NTU standards

- SKU #85522 – Primary Kit
- SKU #85395 – Secondary Kit
- SKU #85548 – Deluxe Kit

FILTER TRACK 660™

#8279

Kit includes 0 and 0.8 NTU standards

- SKU #86341 – Liter Kit
- SKU #86342 – Gallon Kit

PORTABLE METERS

2100 P

#8216

Kit includes 0.0, 20, 100, 800 NTU standards

- SKU #85526 – Primary Kit
- SKU #85407 – Secondary Kit
- SKU #85552 – Deluxe Kit

2100 Q

#8542 White Light

Kit includes 0.0, 20, 100, 800 NTU standards

- SKU #85712 – Primary Kit
- SKU #87010 – Secondary Kit
- SKU #85713 – Deluxe Kit

#8543 Infrared

Kit includes 0.0, 20, 100, 800 FNU standards

- SKU #85714 – Primary Kit
- SKU #85716 – Deluxe Kit

ONLINE METERS

1720 C/D/E

#8247 Full Range

Kit includes 0.0 and 20 NTU standards

- SKU #85558 – Liter Kit
- SKU #85559 – Gallon Kit

#8271 Low Range

Kit includes 0.0 and 1.0 NTU standards

- SKU #86340 – Liter Kit
- SKU #86307 – Gallon Kit

PRODUCTS – PROVEN & TIME TESTED

- Non-toxic; safe handling during use
- Traceable to NIST Primary Standards
- Lot to lot accuracy of 1%
- Easy-to-use: No mixing or dilutions
- Shelf life guaranteed a year from DOS*
- Custom values available upon request



*DOS: Date of Shipment



Individual Standards for Hach Instrumentation

We offer Individual Primary Standards in : 125mL, 250mL, 500mL, 1L, and 1 Gallon.
We also Offer: Sealed Secondary Standards.

2100N/AN (TL2300/2350)	
Parent Bulk Item #	NTU
8215	0, 20, 200, 1000, 4000
8000	0
8098	0.1
8127	0.2
8129	0.3
8290	0.4
8106	0.5
8291	0.6
8292	0.8
8067	1.0
8099	2.0
8687	3.0
8767	4.0
8419	4.5
8100	5.0
8694	6.0
8420	7.0
8768	8.0
8650	10
8422	15
8689	18
8107	20
8688	30
8513	40
8802	50
8691	100
8108	200
8109	1000
8110	4000

2100N IS/AN IS (TL2310/2360)	
Parent Bulk Item #	NTU
8195	0, 20, 200, 1000, 4000
8000	0
8107	20
8508	50
8147	100
8108	200
8509	500
8148	1000
8032	4000

2100P	
Parent Bulk Item #	NTU
8216	0, 20, 100, 800
8000	0
8654	0.1
8293	0.2
8234	0.3
8294	0.4
8115	0.5
8295	0.6
8296	0.8
8461	1.0
8655	2.0
8124	5.0
8111	10
8112	20
8656	40
8113	100
8649	200
8114	800

2100Q WhiteLight	
Parent Bulk Item #	NTU
8542	0, 20, 100, 800
8000	0
8764	0.5
8690	1.0
8765	2.0
8780	4.0
8766	5.0
8534	10
8535	20
8781	40
8536	100
8537	800

2100Q IR	
Parent Bulk Item #	FNU
8543	0, 20, 100, 800
8000	0
8538	10
8539	20
8844	40
8540	100
8541	800

1720 C/D/E	
Parent Bulk Item #	NTU
8271	0, 1.0
8247	0, 20
8000	0
8007	1.0
8070	5.0
8050	20

TU-5 Series	
Parent Bulk Item #	NTU
8000	0
8819	0.1
8820	0.3
8735	0.5
8736	1.0
8821	2.0
8737	10
8738	20
8739	600

2100A	
Parent Bulk Item #	NTU
8201	0.5, 5.0, 40
8202	0.1, 0.5, 5.0

*Please inquire for individual standards

Ratio Instruments	
Parent Bulk Item #	NTU
8212	1.8, 18, 180
8214	1.8, 18, 180, 1800
8000	0
8035	1.8
8037	18
8039	180
8041	1800

GFS CHEMICALS ALSO OFFERS:

- pH Buffers
- Conductivity Standards
- AA/ICP Standards
- UV/Vis Spectrophotometer Standards
- High Purity Acids
- High Purity Solvents

Updated 1/12/24