THERE'S A SAFER WAY TO DETERMINE WATER IN CRUDE OIL.



WATERMARK® Crude Oil Reagent Kits include double ended ampoules containing the KF generator solution - **THE ULTIMATE** in safety & convenience. Packed under vacuum, the generator solution doesn't come out of the ampoule until the second end is snipped.

gfschemicals.com/karl-fischer-reagents

- Crack and Pour for ease and speed no measuring, no glassware to wash.
- Snip Tool reduces glass cuts. Snipped tips remain in tool.
- **Specially formulated** for crude oil applications to ensure complete sample dissolution.

Compatible with all titrator brands.

Watch how it works!



COULOMETRIC REAGENT KITS FOR OILS

SKU: 19026 **ITEM:** 5605

KIT CONTAINS:

8x100 mL bottles of anolyte

• 8x5 mL ampoules of catholyte

One snip tool

Suitable for ASTM D4928





FOR MORE
INFORMATION
service@gfschemicals.com
or call 800.858.9682